FLOW-DIRECTED THERMODILUTION FIBEROPTIC CONTINUOUS CARDIAC OUTPUT PULMONARY ARTERY CATHETER

LIST 011-52509 (Heparinized), 011-52510 (Heparinized), 011-52511

Do not reuse. Reuse negatively impacts performance / sterility potentially resulting in product failure / contamination.

INDICATIONS

For the assessment of the hemodynamic status of a patient, including but not restricted to the following:

- Hemodynamic Pressures
- Cardiac Output
- Oxyhemoglobin Saturation
- Venous Blood Sampling

A secondary indication is for the therapeutic infusion of solutions.

Kinking, excessive manipulation, or grasping with forceps or hemostats can permanently damage the fiberoptics, reducing the transmission of light to and from the blood.

CONTRAINDICATIONS

Use of heparinized catheters is contraindicated in patients with known sensitivity to heparin. Use of any pulmonary artery catheter may be contraindicated in patients with recurrent sepsis, or with a hypercoagulable state where the catheter can serve as a focal point for septic or bland thrombus formation. The Opti-Q™ SvO₂/CCO catheter should not be used in patients weighing less than 20 Kg. In patients this small, the rate of energy delivered via the catheter thermal coil could represent a significant fraction of the patient's total metabolic rate.

SYSTEM OVERVIEW

The Opti-Q™ SvO₂/CCO catheter is intended for use as an integral part of the Oximetrix Oximetry system and the ICU Medical Opti-Q CCO system. The system provides *in vivo* monitoring of oxyhemoglobin saturation (SO₂) without blood sampling and continuous monitoring of cardiac output.

The combined system has four principal components:

1. The Opti-Q SvO₂/CCO catheter is an 8 French catheter with 110 cm of useable length. It contains: the fiberoptics for light transmission; PA, CVP, and distal thermal coil positioning port lumens for pressure readings, sampling or infusion; a thermistor for measurement of cardiac output; a balloon with a volume capacity of 1.5 cc and a thermal coil for continuous measurement of cardiac output.

Markings on the catheter provide an indication of the depth of insertion. A single narrow black line indicates 10 cm and a wide black line indicates 50 cm.

Each catheter is packaged with a disposable black calibration reference assembly for completing automatic calibration prior to insertion. This feature results in accurate SO₂ monitoring immediately upon catheter insertion. *In vivo* calibration is typically not necessary unless the catheter has been inserted into the patient without standardization.

If *in vivo* calibration is necessary, see the appropriate operating manual for proper procedure.

2. The OPTICAL MODULE contains three solid-state light sources (LEDs) and a solid-state light detector. Light from the LEDs is transmitted through the Opti-Q SvO₂/CCO catheter fiberoptics to illuminate the blood at the tip of the catheter. A portion of the light is reflected from the blood and returned through the Opti-Q SvO₂/CCO catheter fiberoptics to the light detector contained within the Optical Module. The reflected light signals are then converted to electrical signals and amplified for transmission to an SO₂ computer.
3. The SO₂ COMPUTER computes values based on the electrical signals transmitted from the Optical Module. These values are continuously displayed in numerical form and are recorded on the unit's display screen trend graph. Visual and audible alarms may signal when the computed oxygen saturation value moves outside the upper and lower limits selected by the operator. The display screen also displays a reflected light intensity signal which provides information on the catheter position and integrity and system function.
4. The CCO COMPUTER, Q-View, controls the catheter thermal coil, detects the resulting blood temperature changes in the pulmonary artery and computes and displays cardiac output. The computer also calculates secondary hemodynamic values.

HEPARIN COATING

The Opti-Q SvO₂/CCO catheter is heparin-coated on all blood contact surfaces.

Although all blood contact surfaces of this catheter are heparin-coated, including inner distal lumen, it is recommended that a constant low rate infusion of heparinized solution be maintained to ensure long term patency of the catheters.

Care should be taken while handling the catheter to ensure that the heparin coating is not removed.

The use of the percutaneous method of insertion is recommended to minimize the possibility of bleeding at the entry site.

CATHETER DESCRIPTION

The OPTI-Q SvO₂/CCO Flow-Directed Thermodilution Fiberoptic Continuous Cardiac Output Pulmonary Artery Catheter is an 8 French, six lumen, radiopaque, polyvinylchloride 110 cm length catheter. The catheter is intended as a monitoring and diagnostic tool providing continuous cardiac output determinations, mixed venous oxygen saturations as well as pressure readings. The catheter is designed to be inserted through an introducer, or via cutdown, into a vein with subsequent positioning of the tip in the pulmonary artery. The catheter is for single use only. See Appendices A & B.

For thermodilution cardiac output only, the OPTI-Q SvO₂/CCO catheter is compatible with other cardiac output computers including Sorenson, Siemens, Marquette Electronics, Nihon Kohden, Mennen, Datascope, and Datex. See Appendix C for thermodilution cardiac output COMPUTATION CONSTANTS.

For other conditions or other cardiac output computers, see the Cardiac Output Computer Operating Manual for the appropriate computation constants, or contact the ICU Medical local sales office.

SETUP

CAUTION: This device should be used by those experienced in the use of flow-directed, balloon-tipped thermodilution catheters.

See the Continuous Cardiac Output System Operating Manual for detailed instructions on use of the system including alarms.

Read the instructions thoroughly before using the catheter. The following procedures assume that the SO₂ computer and optical module are connected and properly stabilized.

NOTE: Inspect sterile packaging first; if damaged, DO NOT USE.

1. Remove outer wrap and leave catheter covered by inner covering.
2. Peel back lift-tab to uncover catheter optical connector. Place optical module in tray recess.
3. Open optical module lid by pulling straight out in direction of arrow, then up. Place catheter optical connector, with the TOP up, in optical module. Close optical module lid.
4. Check black calibration reference assembly. Verify that catheter tip is within reference.
5. Calibrate the SO₂ computer according to manufacturer's instructions.
For Oximetrix 3:
 - a) Press [CAL] key to display the calibration options. Select P for PRE-INSERTION. The following message will appear on the screen as a reminder:
A: CONNECT OPTICS
B: CHECK REFERENCE
 - b) Select Y to initiate pre-insertion calibration. CAL is displayed on the screen and the computer completes calibration procedure within one minute. After successful calibration, the screen displays CAL OK.
 - c) If CAL FAIL is displayed on the screen, repeat the pre-insertion procedure (steps 3- 5). If a second effort fails to produce the CAL OK message, verify that the connection was properly made (step 3) and that the black calibration reference assembly was properly engaged (step 4). Refer to the Oximetry System Operating Manual if still unsuccessful.
6. The oximetry system is now ready for use.
7. Following instructions in the CCO Operating Manual, connect thermal coil connector to thermal coil cable. Setup CCO instrument following instructions in CCO Operating Manual.
8. Pull off remaining inner catheter covering.
9. Using sterile technique, grasp catheter near the entrance of the black calibration reference assembly at the recessed area of the tray and gently lift catheter up and out as shown in Appendix E. Care should be taken in removing the catheter, as the fiberoptics or the balloon may be damaged if the catheter is withdrawn improperly.

10. Test the thermistor for proper function by connecting the thermistor connector to the cardiac output cable and look for an appropriate temperature reading. If there appears to be a thermistor fault, check all connections. If the fault remains, check the cable using the attached test plug. If the fault remains, use a new catheter.
11. Lubricate and test the balloon prior to insertion into a contamination shield or introducer by immersing the balloon in sterile, intravenous solution (e.g., 0.9% Sodium Chloride Injection) and inflating to 1.5 cc while watching for air bubbles.
12. If a contamination shield is used, retest the balloon for proper function after insertion through the shield.
13. Attach sterile stopcocks at the luer hubs of the distal, proximal and distal thermal coil positioning port lumens. Fill the open lumens with an appropriate sterile fluid in order to remove all air, and confirm fluid path and patency of each lumen. Insert the catheter using standard technique and guidelines as follows.

INSERTION

1. The OPTI-Q SvO₂/CCO catheter may be inserted either by cutdown or percutaneously. Routine depth of insertion in normal size adults is 40-50 cm, if inserted via the right jugular vein.
2. Connect the distal lumen and the distal thermal coil positioning port lumen to a pressure monitoring line. For heparinized catheters, a constant low-rate infusion of heparinized solution is essential to reduce the possibility of fibrin deposits over the fiberoptics, to maintain lumen patency, and to provide continuous pressure monitoring. For non-heparinized catheters, use a constant low-rate infusion of saline solution.
3. The catheter is advanced to the right atrium while continuously monitoring pressure at the distal lumen which indicates the tip position.
4. Once in the right atrium, inflate the balloon with 1.5 cc of CO₂ or air and advance the catheter through the right ventricle to the pulmonary artery.
5. A normal pattern of pressure responses while advancing the catheter through the heart chambers should appear as shown in Appendix D.
6. If the attempt to position the catheter fails, deflate the balloon by opening the stopcock and removing the inflation syringe. Pull the catheter back to the right atrium. Repeat Step 4.
7. When the catheter tip is positioned in flowing blood, the system will immediately provide SvO₂ readings.
8. Verify that the thermal coil section extends beyond the tip of the introducer. This assures accurate CCO measurements. The minimal depth of insertion can be determined by adding 29 cm to the length of the introducer.
9. Verify the position of the proximal port. The proximal port should be in the right atrium for optimal thermodilution cardiac output determinations.
10. Verify the position of the distal thermal coil positioning port. The distal thermal coil positioning port should be in the right ventricle for optimal continuous cardiac output determination. (See RV pressure tracing.)
11. Verify the position of the catheter tip and secure the catheter to the patient. The optical module should be secured to the patient or in close proximity to avoid strain or tension on the catheter (see Precautions). **Perform a light intensity calibration** (see Operating Manual).

CAUTION: The catheter should never remain in the wedged position after a reading is taken. To ensure balloon deflation, open the stopcock and remove the inflation syringe. After deflation, reconnect the syringe to prevent inadvertent fluid injection.

PRECAUTIONS

Use during cardiopulmonary bypass

Before initiation of cardiopulmonary bypass, the thermal coil should be deactivated. Consult the Continuous Cardiac Output System Operating Manual for detailed instructions.

Electrocardiographic monitoring during catheter passage is essential, especially where the following conditions are present:

- Complete left bundle branch block where the risk of complete heart block is increased.
- Wolf-Parkinson-White syndrome and Ebstein's malformation where the risk of tachyarrhythmia is present.

NOTE: These risks are less than with standard rigid catheters and in these situations, balloon flotation catheters are specifically indicated to lessen such risks when catheterization is required.

- The following complications have been associated with the use of flow-directed thermodilution catheters: arrhythmias, endocarditis, pneumothorax, arteriovenous fistula formation, damage to vessel intima, right ventricular perforation, pulmonary artery perforation, pulmonary infarction, air embolism and thromboembolic complications.
- Care must be taken not to damage the balloon when advancing the catheter through hemostasis valves of side port catheter/sheath adaptors and rubber seals of contamination shields. Best results are obtained by dilating the valves and rubber seals and moistening the surfaces of the balloon and the seals or valves with a sterile, injectable fluid (for example, normal saline).
- The OPTI-Q SvO₂/CCO catheter is a precision optical instrument and must be handled with care. Kinking, excessive manipulation, or grasping with forceps or hemostats can permanently damage the fiberoptics, reducing the transmission of light to and from the blood. This may result in a low intensity reading, an INTENSITY alarm, and inaccurate SO₂ readings. Kinking or excessive manipulation of the thermal coil section may cause damage and loss of function.
- Patients should not undergo MRI examination with this catheter *in situ*.
- The infusion rate through the distal lumen should not exceed 6ml/hr during CCO monitoring. Furthermore, the combined infusion rate through the proximal and distal thermal coil positioning ports should not exceed 500 ml/hr during CCO monitoring. Higher flow rates may result in erroneous CCO readings.
- If the optical module is not secured directly or in close proximity to the patient, excessive tension will be placed on the catheter, resulting in its movement from the optimal placement. The optical module should be positioned to avoid contact with liquids because fluid entering the catheter-optical module connection may impair light transmission.
- To avoid fibrin formation over the fiberoptics, and to maintain patency of the open lumens, a constant low-rate infusion of heparinized solution is essential. For non-heparinized catheters, use a constant low-rate infusion of saline solution.
- The catheter tip pressure tracings should be monitored continuously. Migration can occur, causing spontaneous wedging which can lead to pulmonary infarction. In the event of spontaneous wedging, reposition the catheter immediately. Chest x-rays should be taken immediately upon insertion and thereafter on a routine basis to ensure that the tip is positioned correctly. Care should also be taken to remove any slack or loops in the catheter as this can also lead to migration.
- The distal thermal coil positioning port tracing should be monitored. The presence of a right ventricular pressure waveform assures that the thermal coil is in the correct place for continuous cardiac output measurements. In the event that a section of the thermal coil migrates out of the right ventricle, turn the power off to the thermal coil and reposition the catheter immediately.
- Knotting can occur with any flow-directed pulmonary artery catheter. Insertion of a guidewire into the PA lumen and manipulation under fluoroscopy may be successful in removing a knot. In some cases, the catheter can be removed if the knot is pulled tight and the catheter is gently and slowly withdrawn. If an introducer is in use, a catheter knot should not be pulled through the introducer. The knot may be withdrawn to the introducer tip and the catheter and introducer removed together.
- Inflation of the balloon to full volume in a small branch of the pulmonary circulation may cause lateral tension sufficient to produce significant damage to the vessel parenchyma, including rupture. Several fatalities resulting from pulmonary artery rupture associated with the use of balloon tipped catheters have been reported.
- The catheter should be well secured to the patient to prevent inward and outward movement. If suture is used, care should be taken that the ties are not too tight as this can damage the fiberoptics or occlude lumens.
- The catheter is supplied sterile and is for single use only. To ensure that the sterility barrier has not been damaged, inspect the outer wrap before using. If there is evidence of damage, do not use the catheter.
- The balloon must not be filled with liquids.
- Although not recommended, if air is used as the inflation medium, the risk of balloon rupture and the possibility of an air embolus entering the arterial system must be assessed. Bacteria-free, filtered CO₂ is the recommended inflation medium if an intracardiac or intrapulmonary shunt is suspected.
- Heparin-coated catheters should not be used on patients with known sensitivity to heparin.
- The Opti-Q SvO₂-CCO catheter is intended for continuous use up to 72 hours.

WARRANTY

ICU Medical warrants that reasonable care has been used in the manufacture of this product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness, since handling, storage, cleaning, or sterilization of this product other than by ICU Medical as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond ICU Medical's control directly affect this product and the results obtained from its use. ICU Medical shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this product other than its replacement. ICU Medical neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this product.

This product or the use thereof is covered by one or more of the following patents: United States: 4,650,327. Canada: 1,086,982; 1,089,252; 1,099,129; 1,277,848. France: 220534. Germany: P3679369.8. Great Britain: 220534. Italy: 220534. Japan: 1940689. Netherlands: 220534.

Store catheters in their individual shipping boxes in a cool area free from exposure to light, dampness, and chemical fumes.

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**MDSS GmbH**
***ICU Medical, Inc.***