



**RE: Use of Transpac® Transducer Devices in MRI**

Dear Clinical Staff;

Thank you for your interest in ICU Medical Transpac Monitoring Sets and Kits. The Transpac transducer is manufactured using tin-coated 24 AWG (American Wire Gage) copper wire in the cable and gold plated copper leads in the transducer chip. None of these materials have ferromagnetic properties which would otherwise cause a magnetic or heat reaction during the MRI procedure. However, these materials do have the potential to cause an artifact and therefore care should be taken to ensure that the transducer and cable are away from the field of view at the time of imaging.

Thank you again for your interest, if there are further questions or concerns please visit our website at [www.icumed.com](http://www.icumed.com), or contact the corporate offices at 949-366-2183.

Technical Services  
ICU Medical, Inc.