

# SENTRY™

SPIKE WITH TEGO™ CONNECTOR




SENTRY™



ICU Medical, Inc.  
♥ CRITICAL CARE SYSTEMS

# SENTRY

SPIKE WITH TEGO™ CONNECTOR




The Sentry spike with Tego connector adds an additional level of safety to your contrast delivery system, creating a microbiologically and mechanically closed system which protects the spike assembly from contamination. This patented connector along with proper disinfecting swabbing procedures will reduce the possibility of patient infection and cross contamination dramatically.



## UNIQUE ADVANTAGES

The name TEGO originates from the Latin word conTEGO which means to protect or defend and that is just what the TEGO connector does.

The Tego connector is removable for ease of use if initial priming of contrast spike is desired.




The user selectable vent on the spike is for use with bottles or bags of contrast.



The simple yet efficient design replaces the use of a standard stopcock.

When used with any contrast management system, it provides the proper length of tubing to give adequate head pressure so that the contrast burette fills and refills easily to the proper level.



PART NO	UNITS PER BOX	UNITS PER CASE
42343-04, AG-7350	25	100